




ACM-ERK410 ICU BED WITH 4 MOTORS



- Backrest, Legrest, Up-Down, Trendelenburg / Anti Trendelenburg positions are adjusted by motors with side rail control panels,
- Cardiac position, Fowler position, Examination position, Auto Regression,
- Head and Foot parts are made of PP and they are easily removable.
- 4 pcs side Rails made of PP,
- 4 pcs 150mm central locking castors,
- Mattress platform is made of cleanable ABS.
- 4 bumpers in every corner.
- All corners IV rod places available, 1 pcs IV rod.
- Angle indicator for backrest
- Drainage bag holders
- Battery (optional)
- Nurse control unit (optional)
- Overall sizes : 2130 x 1030 x (430-780) mm
- Mattress platform : 1950 x 850 mm



4 MOTORS

CPR

LOCKED WHEELS

KG 230