



ACM-ES500 EMERGENCY STRETCHER



- Head part of the stretcher can move up to 62 degrees by two pieces gas springs,
- Height, Trendelenburg/reverse Trendelenburg are adjusted by hydraulic pumps,
- Laying surface is made of x-ray permeable HPL material
- 200 mm central brake castors,
- 125 mm 5th castor,
- X-ray cassette holder full length of the stretcher,
- Protective rubber pumbers on 4 corners,
- Sizes: 2170 mm x 840 mm x (560-870 mm),
- Safe loading capacity 300 kg.,

