



ACM-GC300 GYNECOLOGICAL CHAIR WITH THREE MOTORS



- Gynaecology Chair – Standard Column is designed for gynecology, examination, urology and minor surgical operations.
- The chair consists of 3 main parts. Back rest, seat part and leg rest.
- The chair makes back up&down, chair up&down and trendelenburg&reverse trendelenburg movements by 3 DC motors with using handset.
- Height can be adjusted by handset.
- It moves up&down vertically.
- Steel construction is painted with electrostatic powder coating and protected with galvanize against corrosion.
- And the construction is covered with anti-bacterial and easy cleanable ABS plastic.
- So that, the chairs seems more aesthetically.
- ABS covering designed to keep itself in the upholstery and it blocks the dust which may go into the upholstery.
- Drain pan does not effected with table movements during the operations, it stays stabile. There are stainless steel rails at the two sides of the table.
- With these rails, all accessories can be moved to all directions easily (up/down/right/left).
- Leg supports are easily removable and surgeons can get the supports manually to all desired positions.

