



ACM-SST003 SCRUB SINK TRIPLE



- Made of 304 grade stainless steel for use in sterile areas
- Flexible production in different designs and dimensions,
- For entire sterilization required places they serve as the scrub station for surgeons and other OR staff to scrub their hands and arms before a surgical or invasive procedure.
- Knee and photocell operated water and soap flow system,
- Stainless steel gooseneck faucet with filter system to leach the particles
- User friendly curved basin design to reduce the splashing of soiled water
- Digital water flow timer, controlled by a programmable electronic timer
- Optional thermostatic mixing valve to adjust water temperature and UV
- Lamp and mirror
- Stainless steel adjustable feet.
- Single, Double and Triple models are available

