



ACM-BLD04 BLOOD DONOR CHAIR WITH 4 MOTORS



- Blood Donor Chair is designed for blood transfusion, haemodialysis, eye operations, skin care, body care, epilation, ozone therapy, oncology and other medical aesthetic operations.
- The chair consists of 3 main parts. Backrest, seat part and legrest.
- The chair makes back up & down, height up & down, legrest up & down and Trendelenburg & reverse Trendelenburg movements by 4 DC motors with using handset.
- Chair is covered with cleanable, waterproof material, also the chair is comfortable and orthopedically designed.