



### ACM-CC001E CRASH CART EPOXY PAINTED



- 55x77x110 cm.,
- Emergency Cart, Crash Cart,
- Electrostatic Powder Coated,
- 5 Drawers with Central Lock System,
- Push Handle,
- Telescopic Drug Preparation Table,
- Height Adjustable IV pole,
- Swiveling Defibrillator Table,
- CPR Board (Plexiglass),
- Wastebin,
- Oxygen Bottle Holder,
- Side Rails for Accessories (2 PCS),
- 1.and 2. Drawers Include 20 Dividers.