



ACM-ERK210 HOSPITAL BED WITH 2 MOTORS AND TRENDELENBURG



- Patient bed with two motors
- Hand remote control unit
- Backrest and legrest adjusted by motor
- Trendelenburg position is adjusted by gas piston
- Removable pp head and foot boards
- PP side rails
- 125mm with two brake castors
- Metal or abs mattress platform (optional)
- All metal parts are epoxy painted with serum holder
- Bumpers on 4 corners of the bed
- Overall sizes: 2130mm x 1030mm x 430mm
- Mattress platform sizes: 1900mm x 850mm



2 MOTORS

KG 230