



ACM-ERK310 HOSPITAL BED WITH 3 MOTORS AND TRENDLENBURG



- Backrest, Legrest, Up-Down, positions are adjusted by motors with Handset,
- Cardiac position, Fowler position, Examination position, Auto Regression,
- Trendelenburg position is adjusted by gas piston
- Head and Foot parts are made of PP and they are easily removable.
- Manual CPR for head part in case of emergency (optional)
- 4 pcs side Rails made of PP,
- 4 pcs 125 mm two with brake castors,
- Mattress platform is made of cleanable ABS plastic.(optional)
- Overall sizes : 2140 x 1040 x (410-700) mm
- 4 bumpers in every corner.
- All corners IV rod places available, 1 pcs IV rod.
- Angle indicator
- Drainage bag holders



3 MOTORS **KG 230**